

# FlexiFit™ 432

FULL FACE MASK

## Specification Sheet


---

A Full Face Delivery Mask




*Our second generation Full Face Mask incorporates a soft Foam Cushion that works together with the highly contoured Silicone Seal to provide improved sealing technology and a more comfortable fit.*

# Full Face Delivery - FlexiFit™ 432




*The FlexiFit™ 432 Full Face Mask offers an improved seal and more comfortable fit with the inclusion of a soft Foam Cushion that interacts with the Silicone Seal. This alternative to nasal or oral delivery comes fully assembled and ready for use.*

## UNIQU MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit™ Technology Auto-contouring fit

*With Foam Cushion* FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments. The Foam Cushion improves the sealing technology and comfort.



### Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



### Under-Chin Design Extra stability

An Under-Chin fit provides stability by providing a constant point of reference and maintains the seal even with movement of the jaw during the night.



### Efficient Bias Diffuser Quiet usage

*Maintenance Free*  
This advanced air diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.




### Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



### Crown Strap Extra stability

Complements the Under-Chin Design to maximize stability and eliminate the need for a chin strap.



### Easy-Release Strap One-step removal

Simply pull the Easy-Release Cord for one-step removal of headgear. No ongoing headgear adjustments required.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- The FlexiFit™ 432 Full Face comes in three separate sizes, Small, Medium and Large
- User friendly Fitting Guide

REORDER CODE: HC432

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE:

**International** P O Box 14 348, Panmure, Auckland 1134, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

**USA** Tel: 1800 446 3908 +1 949 470 3900 Fax: +1 949 470 3933 **Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22  
**China** Tel: +86 20 3205 3486 Fax: +86 20 3205 3495 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99  
**India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205  
**UK (EU Authorised Representative)** Tel: +44 1628 626 136 Fax: +44 1628 626 146 **Nthn Ireland** Tel: 0800 132 189  
**Netherlands** Tel: +31 78 644 0924 Fax: +31 78 644 0914 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Switzerland** Tel: 0800 83 47 63  
**USA** Fax: +81 3 3661 7206